## HEINEn suorat oftalmoskoopit

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| Guidance and manufacturer's declaration – electromagnetic emissions                                                                                                          |                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment. |                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Emission test                                                                                                                                                                | Compliance     | Electromagnetic environment – Guidelines                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| RF emissions CISPR11                                                                                                                                                         | Group 1        | The device uses RF energy only for its internal function. Therefore, RF-emission is very low and it is unlikely that any interference in nearby electronic equipment.                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| RF emissions CISPR 11                                                                                                                                                        | Class B        | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. <b>Warning:</b> This device is intended only for use by medical professionals. This is a device of class A CISPR 11 in the domestic environment, this device may cause radio interference, so that it may be necessary in this case, to take appropriate remedial measures, as e.g. orientation, new arrangement or shielding of the device or restrict the connection to the site. |
| Harmonic Emissions IEC 61000-3-2                                                                                                                                             | not applicable |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3                                                                                                                        | Passed         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |

| Guidance and manufacturer declaration - Electromagnetic immunity                                                                                                             |                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                                 |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment. |                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                                 |
| Immunity test                                                                                                                                                                | IEC 60601 test level                                                                                                                                                                                                                                                         | Compliance level                                                                                                                                                                                                                                                             | Electromagnetic environment – Guideline                                                                                                                                                                                                                                         |
| Electrostatic discharge (ESD) acc. to IEC 61000-4-2                                                                                                                          | ± 6 kV contact discharge<br>± 8 kV air discharge                                                                                                                                                                                                                             | ± 6 kV contact discharge<br>± 8 kV air discharge                                                                                                                                                                                                                             | Floors should be wood, concrete covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.                                                                                                                       |
| Electrical fast transient/ burst IEC 61000-4-4                                                                                                                               | ± 2 kV for mains cables<br>± 1 kV for input and output lines                                                                                                                                                                                                                 | ± 2 kV for mains cables<br>± 1 kV for input and output lines                                                                                                                                                                                                                 | The supply voltage quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                      |
| Surge IEC 61000-4-5                                                                                                                                                          | ± 1 kV voltage phase – phase,<br>± 2 kV voltage phase – earth                                                                                                                                                                                                                | ± 1 kV voltage phase – phase<br>± 2 kV voltage phase – earth                                                                                                                                                                                                                 | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                             |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11                                                                          | <5% U <sub>T</sub> , (>95% dip in U <sub>T</sub> ) for 1/2 period<br>40% U <sub>T</sub> , (60% dip in U <sub>T</sub> ) for 5 periods<br>70% U <sub>T</sub> , (30% dip in U <sub>T</sub> ) for 25 periods<br><5% U <sub>T</sub> , (>95% dip in U <sub>T</sub> ) for 5 seconds | <5% U <sub>T</sub> , (>95% dip in U <sub>T</sub> ) for 1/2 period<br>40% U <sub>T</sub> , (60% dip in U <sub>T</sub> ) for 5 periods<br>70% U <sub>T</sub> , (30% dip in U <sub>T</sub> ) for 25 periods<br><5% U <sub>T</sub> , (>95% dip in U <sub>T</sub> ) for 5 seconds | Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered by a UPS (uninterruptible power supply) or a battery. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8                                                                                                                      | 3 A/m                                                                                                                                                                                                                                                                        | 3 A/m                                                                                                                                                                                                                                                                        | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                       |
| Comment: U <sub>T</sub> is the a.c. supply voltage prior to application of the test level.                                                                                   |                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                                 |